# PRONTO® .035"

Extraction Catheter

## PRONTO® .035" Extraction Catheter Model 5040

### Instructions For Use

#### USA CAUTION

Federal (USA) law restricts this device to sale by or on the order of a physician.

#### CAUTION

The PRONTO .035" Extraction Catheter procedure should be performed by physicians with adequate training in the use of the device.

#### DEVICE DESCRIPTION

The PRONTO .035" Extraction Catheter (PRONTO .035") is a dual lumen, over-the-wire (OTW) catheter with related accessories. The catheter is designed to be delivered through a 10F or larger introducer sheath over a 0.035" guidewire. The larger lumen allows for the removal of thrombus by use of the included syringe through the extension line and stopcock. The catheter has a rounded distal tip with a protected extraction lumen to facilitate advancement of the catheter into the blood vessel and to maximize extraction of thrombus through the extraction lumen.

The catheter has a radiopaque marker band located approximately 4mm from the distal tip. The proximal end of the catheter incorporates a hemostatic Y-junction that allows for the attachment of the catheter to the included extension line, stopcock and syringe and can be tightened onto the guidewire to prevent blood leakage.

When using the PRONTO .035" for fluid delivery, the maximum flow rate recommended for the device is as follows (see Table 1):

**TABLE 1: Pronto .035" Extraction Catheter Flow Rate**

Maximum Recommended Catheter Flow Rate	
Saline	1.00ml/s
76% Ionic Contrast Media	1.00ml/s

#### INDICATIONS

The PRONTO .035" is indicated for:

- the removal/aspiration of embolic material (thrombus/debris) from vessels of the arterial system
- the removal/aspiration of embolic material (thrombus/debris) from vessels of the deep venous system
- to infuse/deliver diagnostic or therapeutic agents

#### CONTRAINDICATIONS

The PRONTO .035" is contraindicated in:

- vessels <4.0mm in diameter
- the removal of fibrous, adherent or calcified material (e.g., chronic clot, atherosclerotic plaque, chronic pulmonary embolism)
- use in the right heart or pulmonary arteries during active cardiopulmonary resuscitation

#### WARNINGS

The PRONTO .035" is supplied sterile and is for single use only. Do not resterilize and/or reuse the device.

Do not use the syringe, extension line or stopcock inside the human body.

If flow into the syringe stops or is restricted, DO NOT discontinue suction; DO NOT attempt to flush the extraction lumen while the catheter is still inside the patient's vasculature. Intravascular thrombus delivery, a thromboembolic event and/or serious injury, or death may result. Remove the catheter

while the suction is on, and once the catheter is outside the patient, either flush the extraction lumen or use a new catheter.

The safety and effectiveness of the Pronto .035" has NOT been established for use as embolic protection or for coronary thrombectomy.

#### COMPLICATIONS

As with all catheterization procedures, complications with the PRONTO .035" may occur. These may include:

- local or systemic infection
- access site bleeding, including local hematoma
- vessel spasm, dissection, perforation or rupture
- venous valve injury
- arrhythmias
- ventricular perforation
- hemoptysis
- pleural effusion
- arterial or venous thrombosis
- distal embolization of blood clots and plaque, pulmonary embolism, or myocardial or pulmonary infarction
- arteriovenous fistula formation
- catheter fracture with tip separation and distal embolization

#### PRECAUTIONS

The PRONTO .035" deployment procedure should be performed by physicians thoroughly trained in percutaneous, intravascular techniques and procedures.

Do not use the PRONTO .035" if the packaging has been damaged.

Inspect the catheter prior to use for any bends or kinks. Do not use a damaged catheter because vessel damage and/or inability to advance or withdraw the catheter may occur.

Check that all fittings are secure so that air is not introduced into the extension line or the syringe during extraction.

Exercise care in handling of the catheter during a procedure to reduce the possibility of accidental breakage, bending or kinking.

When the catheter is in the body, it should be manipulated only under fluoroscopy. Do not attempt to move the catheter without observing the resultant tip response.

Never advance or withdraw an intravascular device against resistance until the cause of the resistance is determined by fluoroscopy. Movement of the catheter or guidewire against resistance may result in separation of the catheter or guidewire tip, damage to the catheter, or vessel perforation.

Precautions to prevent or reduce clotting should be taken when any catheter is used in the vascular system. Use of systemic heparinization and heparinized sterile solutions should be considered.

#### CLINICAL PROCEDURE

The following instructions provide technical direction but do not obviate the necessity of formal training in the use of the PRONTO .035". The techniques and procedures described do not represent ALL medically acceptable protocols, nor are they intended as a substitute for the physician's experience and judgment in treating any specific patient.

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Each PRONTO .035" includes the following components:

- Single-use, disposable catheter
- Extension line with attached roller clamp and 3-way stopcock
- 60ml syringe with locking plunger

Other materials required but not provided are:

- 10F or greater introducer sheath
- Sterile heparinized saline (for system flushing)

#### PREPARATIONS FOR USE

1. sterile technique, open the pouch and transfer the packaging card with catheter and accessories into the sterile field.
2. Remove the catheter and accessories from the packaging card and inspect the catheter and extension tubing for any bends or kinks.
3. Remove the packaging stylet and thoroughly back-flush the guidewire lumen with heparinized saline.
4. Draw 10ml of heparinized saline into the 60ml syringe and ensure that all connections between the syringe, stopcock and extension line are secure and that the yellow cap is in place on the stopcock side port. Connect the extension line to the side port of the tuohy adapter on the catheter hub.
5. Flush all 10mls of heparinized saline through the catheter and turn the stopcock "off" to the syringe. Fully retract the syringe plunger and twist to lock in place.
6. Set the roller clamp so that the roller is in the "off" position (toward the syringe).  
  
Precaution: Check that all fittings are secure so that air is not introduced into the extension line or syringe during extraction.

#### DEPLOYMENT PROCEDURE

Deploy the PRONTO .035" according to the following steps:

1. Cannulate the vessel using standard technique and a properly prepared 10F or greater introducer sheath.
2. Insert the guidewire into the desired vascular space.
3. Back-load the PRONTO .035" onto the guidewire and advance it over the guidewire until the wire exits the proximal end of the Y-junction.
4. Secure the proximal segment of the guidewire and insert the PRONTO .035" through the introducer sheath and into the desired vascular space. Stop advancement of the catheter if any resistance is encountered.

**WARNING: Never advance or withdraw an intravascular device against resistance until the cause of the resistance is determined by fluoroscopy. Movement of the catheter or guidewire against resistance may result in separation of the catheter or guidewire tip, damage to the catheter, or vessel perforation.**

5. After fluoroscopically confirming catheter position, open the stopcock so it is open to the extension line.
6. Open the roller clamp until blood enters the syringe at the desired rate.
7. Slowly advance the catheter distally into the thrombus. Blood will enter the syringe until the vacuum is gone. Should aspiration not begin filling the syringe within 5 seconds, remove the catheter while still extracting. Outside the patient, either flush the extraction lumen or use a new catheter.

**WARNING: If flow into the syringe stops or is restricted, do NOT attempt to flush the extraction lumen while the catheter is still inside the patient's vasculature. Intravascular thrombus delivery, a thromboembolic event and/or serious injury, or death may result. Remove the catheter while the suction is on, and once the catheter is outside the patient, either flush the extraction lumen or use a new catheter.**

8. After completing the extraction process, turn the stopcock so that it is "off" to the catheter and remove the catheter.

#### PACKAGING & STORAGE

The PRONTO .035" catheter has been sterilized with ethylene oxide.

#### LIMITED WARRANTY

Vascular Solutions, Inc. warrants that the Pronto .035" Extraction Catheter is free from defects in workmanship and materials prior to the stated expiration date. Liability under this warranty is limited to refund or replacement of any product which has been found by Vascular Solutions, Inc. to be defective in workmanship or materials. Vascular Solutions, Inc. shall not be liable for any incidental, special, or consequential damages arising from the use of the Pronto .035" Extraction Catheter. Damage to the product through misuse, alteration, improper storage or improper handling shall void this limited warranty.

No employee, agent or distributor of Vascular Solutions, Inc. has any authority to alter or amend this limited warranty in any respect. Any purported alteration or amendment shall not be enforceable against Vascular Solutions, Inc.

THIS WARRANTY IS EXPRESSLY IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE OR ANY OTHER OBLIGATION OF VASCULAR SOLUTIONS, INC.

#### PATENTS AND TRADEMARKS

International and U.S. patents pending.

Pronto® is a registered trademark of Vascular Solutions, Inc.

See the International Symbols Glossary on page 8.

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PRONTO® es una marca registrada de Vascular Solutions, Inc.