TWIN-PASS® Dual Access Catheters Instructions For Use

USA CAUTION

Federal (USA) law restricts this device to sale by or on the order of a physician.

CAUTION

The TWIN-PASS dual access catheters should be used by physicians with adequate training in the use of the device. The device is provided sterile for single use only. Reuse of single-use device creates a potential risk of patient or user infections. Contamination of the device may lead to illness or serious patient injury.

DEVICE DESCRIPTION

The TWIN-PASS dual access catheters are dual lumen catheters designed for use in the arterial vasculature. They consist of a rapid exchange delivery lumen on the distal segment, and an OTW lumen that runs the length of the catheter. The catheters are available in different lumen diameter configurations. The distal tip of each catheter is identified by a radiopaque marker band, and some configurations contain a second proximal marker band identifying the end of the OTW lumen. The TWIN-PASS catheters are provided with a stiffening mandrel that, when loaded into the OTW lumen, provides support and pushability during catheter insertion.

Product Features

Model Number	5200	5230
Max Guidewire RX lumen, in / mm	0.014" / 0.36mm	0.014" / 0.36mm
Max Guidewire OTW lumen, in / mm	0.014" / 0.36mm	0.018" / 0.46mm
RX lumen I.D., in / mm	0.017" / 0.43mm	0.017" / 0.43mm
OTW lumen I.D., in / mm	0.017" / 0.43mm	0.023" / 0.58mm
Proximal O.D., F	3F	3.5F
Distal O.D., F	3F	3.5F
Tip O.D., F	1.9F	1.9F
Minimum Guide Catheter I.D., in / mm	0.058" / 1.47mm	0.066" / 1.68mm

Max Recommended Flow Rate

Model Number	5200	5230
Saline	0.31ml/sec	1.24ml/sec
76% Ionic Contrast	0.047ml/sec	0.163ml/sec

INDICATIONS

The TWIN-PASS catheters are intended to be used in conjunction with steerable guidewires in order to access discrete regions of the coronary and peripheral arterial vasculature, to facilitate placement and exchange of guidewires and other interventional devices, for use during two guidewire procedures and to subselectively infuse/deliver diagnostic or therapeutic agents.

The TWIN-PASS .023" catheter (Model 5230) is also intended to be used to measure intra-arterial pressure within the peripheral and coronary vasculature.

CONTRAINDICATIONS

The Twin-Pass catheter is contraindicated for high pressure injections and for use in the cerebral vasculature.

WARNINGS

Do not advance the TWIN-PASS catheter without a guidewire in place through the RX lumen. Advancement of the catheter without a guidewire in the RX lumen may result in intimal damage, arterial dissection, or perforation.

Do not reuse, reshape or re-sterilize the device. Re-sterilization or reshaping may change the physical characteristics of the material and should not be attempted.

Do not exceed the maximum recommended flow rate when injecting through the TWIN-PASS catheter. Catheter rupture and arterial injury could result.

If a guidewire of appropriate diameter cannot be passed through the catheter, do not attempt to resolve the blockage by flushing the catheter in vivo. Catheter rupture and arterial injury could result. Identify and resolve the cause of the blockage, or replace the catheter with a new one.

Never advance or withdraw an intravascular device against resistance until the cause of the resistance is determined by fluoroscopy. Movement of the catheter or guidewire against resistance may result in separation of the catheter

or guidewire tip, damage to the catheter, or vessel perforation.

COMPLICATIONS

As with all catheterization procedures, complications may occur when using the TWIN-PASS catheter. These may include:

- local or systemic infection
- intimal disruption
- arterial dissection
- perforation and vessel rupture
- arterial thrombosis
- distal embolization of blood clots and plaque
- myocardial infarction
- arterial spasm
- catheter fracture with tip separation and distal embolization

PRECAUTIONS

The TWIN-PASS catheter deployment procedure should be performed by physicians thoroughly trained in percutaneous, intravascular techniques and procedures.

Do not use the TWIN-PASS catheter if the packaging has been damaged.

Inspect the catheter prior to use for any bends or kinks. Do not use a damaged catheter because vessel damage and/or inability to advance or withdraw the catheter may occur.

Both catheter lumens must be flushed with sterile, heparinized saline prior to use.

Precautions to prevent or reduce clotting should be taken when any catheter is used in the vascular system. Use of systemic heparinization and heparinized sterile solution should be considered.

Exercise care while handling the catheter during a procedure to reduce the possibility of accidental breakage, bending or kinking.

Excessive tightening of a hemostatic valve onto the catheter shaft may result in damage to the guidewire lumen, difficulty while inserting the catheter or guidewires, or damping of measured pressures.

When the catheter is in the body, it should be manipulated only under fluoroscopy. Do not attempt to move the catheter without observing the resultant tip response.

Check that all fittings and attachments are secure and that air is not introduced into the system. The presence of an air bubble or the use of an inappropriate extension line (excessive length, ID or compliance) can cause damping of the pressure signal.

CLINICAL PROCEDURE

The following instructions provide technical direction but do not obviate the necessity of formal training in the use of the TWIN-PASS catheter. The techniques and procedures described do not represent ALL medically acceptable protocols, nor are they intended to be a substitute for the physician's experience and judgment in treating any specific patient.

When using the TWIN-PASS catheter as a diagnostic aid, all available data, including the patient's signs and symptoms and other diagnostic test results, should be considered before determining a specific treatment plan.

Each TWIN-PASS dual access catheter includes the following components:

- Single-use disposable catheter
- Stiffening mandrel
- Dispenser coil with flushing luer

Other materials required but not provided are:

- Guiding catheter with an I.D. of at least 0.058" / 1.47mm fitted with a rotating hemostatic valve (RHV) (Tuohy-Borst type)
- Appropriate guidewires
- 10ml syringe (for flushing the dispenser coil and the catheter lumen)
- Sterile heparinized saline (for system flushing)
- Transducer (if measuring pressure)

PREPARATIONS FOR USE

1. Carefully inspect the TWIN-PASS catheter packaging and components for damage prior to use. Utilizing sterile technique, remove the TWIN-PASS catheter

English/Instructions for Use	1
Český/Návod k použití	2
Dansk/ Brugsanvisning	3
Nederlands/Gebruiksaanwijzing.....	5
Suomi/Käyttöohjeet	6
Français/Mode d'emploi	7
Deutsch/Gebrauchsanweisung	9
Ελληνικά/Οδηγίες χρήσης	10
Italiano/Istruzioni per l'uso.....	11
Norsk/Bruksanvisning	13
Polski/Instrukcja stosowania	14
Português/Instruções de Utilização.....	15
Español/Instrucciones de uso	16
Svenska/Bruksanvisning	18
Türkçe/Kullanım Talimatları	19

Vascular Solutions, Inc.
6464 Sycamore Court
Minneapolis, MN 55369 USA
(888) 240-6001 USA
(763) 656-4300
(763) 656-4250
www.vasc.com

Vascular Solutions Zerusa Limited
208 Business Innovation Centre
NUI Galway
Newcastle Road
Galway Ireland
(+353) 91 861611
(+353) 91 861612

- dispenser coil from its packaging and transfer it to the sterile field.
- Remove the stiffening mandrel from the dispensing coil. DO NOT DISCARD.
 - Attach a 10ml syringe filled with sterile heparinized saline to the luer-lock guidewire entry port of the TWIN-PASS catheter and thoroughly flush the catheter.
 - Completely flush the coil to activate the hydrophilic coating on the TWIN-PASS catheter (if applicable).
 - Insert the stiffening mandrel through the luer-lock and into the TWIN-PASS catheter and lock it in place.
 - Remove the TWIN-PASS catheter from the dispensing coil and inspect it for any bends or kinks.
 - Remove the packaging mandrel from the rapid exchange lumen of the TWIN-PASS catheter while under sterile saline.

DEPLOYMENT PROCEDURE

The following TWIN-PASS catheter deployment steps assume a standard PTCA protocol using the following items: a guiding catheter, an inserted guidewire, an exchange length guidewire of the appropriate diameter to be delivered through the OTW lumen.

As with any interventional procedure, proper anticoagulation and anti-platelet therapy should be administered prior to beginning.

Note: Familiarity with traditional long and short guidewire exchange techniques is required for successful deployment of the TWIN-PASS catheter and the delivery of a second guidewire.

TWIN-PASS DEPLOYMENT STEPS

- Backload the rapid exchange segment of the TWIN-PASS catheter onto the proximal end of the guidewire that is already in place in the distal vasculature.

WARNING: Do not advance the TWIN-PASS catheter without a guidewire in place through the RX lumen. Advancement of the catheter without a guidewire in the RX lumen may result in intimal damage, arterial dissection, or perforation.

- Carefully advance the catheter into the desired distal vascular space.

WARNING: Never advance or withdraw an intravascular device against resistance until the cause of the resistance is determined by fluoroscopy. Movement of the catheter or guidewire against resistance may result in separation of the catheter or guidewire tip, damage to the catheter, or vessel perforation.

- Slowly remove the stiffening mandrel.
- To deliver an exchange-length guidewire through the OTW lumen:

- Load the guidewire into the luer-lock of the TWIN-PASS catheter.
- Advance the guidewire until it exits the OTW lumen into the distal vascular space.
- Fix both guidewires using standard guidewire exchange techniques and carefully withdraw the TWIN-PASS catheter until the distal tip exits the hemostatic valve and both wires can be secured.

- To deliver fluid through the OTW lumen:
 - Create an air-free, fluid-to-fluid connection between the OTW lumen and the injection syringe.
 - Aspirate prior to injection to remove any air bubbles remaining in the catheter.

Warning: Do not exceed the maximum recommended flow rate when injecting through the TWIN-PASS catheter. Catheter rupture and arterial injury could result.

- To measure intra-arterial pressure through the OTW lumen (TWIN-PASS 0.023" only):
 - Connect the OTW lumen to a pressure transducer using fluid-to-fluid connection.
 - Remove all air bubbles from the catheter, extension tubing (if used) and transducer.

Precaution: Check that all fittings and attachments are secure and that air is not introduced into the system. The presence of an air bubble or the use of

an inappropriate extension line (excessive length, ID, or compliance) can cause damping of the pressure signal.

When not in use during this procedure, wipe the TWIN-PASS catheter with a sterile gauze pad saturated with heparinized saline, flush the OTW lumen well, reload the stiffening mandrel, and store in the dispensing tube in a saline bath.

PACKAGING & STORAGE

The TWIN-PASS catheter has been sterilized with ethylene oxide.

STERILE EO

Handle with care.

LIMITED WARRANTY

Vascular Solutions, Inc. warrants that the TWIN-PASS catheter is free from defects in workmanship and materials prior to the stated expiration date. Liability under this warranty is limited to refund or replacement of any product that has been found by Vascular Solutions, Inc. to be defective in workmanship or materials. Vascular Solutions, Inc. shall not be liable for any incidental, special, or consequential damages arising from the use of the TWIN-PASS catheter. Damage to the product through misuse, alteration, improper storage, or improper handling shall void this limited warranty.

No employee, agent, or distributor of Vascular Solutions, Inc. has any authority to alter or amend this limited warranty in any respect. Any purported alteration or amendment shall not be enforceable against Vascular Solutions, Inc.

THIS WARRANTY IS EXPRESSLY IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE OR ANY OTHER OBLIGATION OF VASCULAR SOLUTIONS, INC.

PATENTS AND TRADEMARKS

International and U.S. patents pending.

Twin-Pass® is a registered trademark of Vascular Solutions, Inc.

See the International Symbols Glossary on page 23.

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Twin-Pass® Vascular Solutions, Inc.'in tescilli bir ticari markasıdır.