Gandras[®] Catheter

Gandras[®] Visci-G Catheter

GANDRAS[®]/GANDRAS Visci-G Catheter Instructions For Use

USA CAUTION

Federal (USA) law restricts this device to sale by or on the order of a physician.

CAUTION

The GANDRAS/GANDRAS Visci-G catheter should be used by physicians with adequate training in angiographic procedures and the device. The GANDRAS/GANDRAS Visci-G is provided sterile for single use only. Reuse of single-use device creates a potential risk of patient or user infections. Contamination of the device may lead to illness or serious patient injury.


DEVICE DESCRIPTION

The GANDRAS/GANDRAS Visci-G catheter is a hydrophilically coated angiographic catheter for use in performing angiographic procedures. Catheters for use in angiographic procedures are available in a variety of French sizes, endhole sizes, lengths, materials and tip configurations. The choice of a specific catheter for an angiographic procedure should be based on the physician's judgment and experience.


The GANDRAS/GANDRAS Visci-G catheter is to be deployed through a standard introducer sheath of appropriate size and will accommodate a standard 0.035" (0.89mm) diameter guidewire.

A single-use, peelable tip straightener is supplied with the GANDRAS/GANDRAS Visci-G catheter.

GANDRAS

GANDRAS Visci-G

INDICATIONS

The GANDRAS/GANDRAS Visci-G catheter is designed to be used for delivering embolic materials and radiopaque media to selected sites in the vascular system. Diagnostic, embolic, or therapeutic agents are to be used in accordance with specifications outlined by the manufacturer.

CONTRAINDICATIONS

The GANDRAS/GANDRAS Visci-G catheter is contraindicated for use in the cerebral vasculature.

WARNINGS

Never advance or withdraw an intravascular device against resistance until the cause of the resistance is determined by fluoroscopy. Movement of the catheter against resistance may

English/Instructions for Use.....	1
Český/Návod k použití	3
Dansk/Bruksanvisning	5
Nederlands/Gebruiksaanwijzing	8
Suomi/käyttöohjeet.....	10
Français/Mode d'emploi	12
Deutsch/Gebrauchsanweisung	14
Ελληνικά/Οδηγίες χρήσης	17
Italiano/Istruzioni per l'uso.....	19
Norsk/Bruksanvisning	21
Polski/Instrukcja użytkowania.....	23
Português/Instruções de Utilização.....	26
Español/Instrucciones de uso	28
Svensk/Bruksanvisning på svenska	30

Vascular Solutions, Inc.
6464 Sycamore Court
Minneapolis, MN 55369
USA
(888) 240-6001 USA
(763) 656-4300
(763) 656-4250
www.vascularsolutions.com

TheraGenesis GmbH
Bruehlstrasse 50
76297 Stutensee
Germany
+49 (0) 724 9952 063
+49 (0) 724 9952 064

result in separation of the catheter tip, damage to the catheter or vessel perforation.

Do not exceed the maximum pressure rating listed on the label of the GANDRAS/GANDRAS VISCI-G catheter for any injection. Exceeding the pressure rating may result in separation of the catheter, embolization, vessel damage or delivery of the therapeutic agent to an unintended site. The maximum recommended pressure rates for the device are as follows (see Table 1):

Table 1: Recommended Pressure Rates

Model	Catheter size	Length	Tip length (A)	Maximum psi rating	Maximum kPa rating
5580	5F	82 cm	2 cm	600	4137
5581	5F	84 cm	4 cm	600	4137
5582	5F	86 cm	6 cm	600	4137
5585	5F	75 cm	1 cm	600	4137

Check that all fittings and attachments are secure so that air is not introduced into the catheter during use.

Do not attempt to resolve a lumen blockage by flushing the catheter *in vivo*. Catheter rupture and arterial injury could result. Identify and resolve the cause of the blockage, or replace the catheter with a new one.

Do not use the GANDRAS/GANDRAS VISCI-G catheter to deliver embolic materials that require a lumen size greater than 0.037" (0.94mm). Using the GANDRAS/GANDRAS VISCI-G catheter in this situation may result in flow obstructions within the catheter and/or catheter damage.

COMPLICATIONS

As with all angiographic procedures, complications with the GANDRAS/GANDRAS VISCI-G catheter may occur. These may include:

- local or systemic infection
- stroke and death
- intimal disruption
- arterial dissection
- perforation and vessel rupture
- arterial or venous thrombosis
- arterial or vessel spasm
- catheter fracture with tip separation and distal embolization
- delivery of embolic material to unintended anatomical location
- arteriovenous fistula formation

PRECAUTIONS

The GANDRAS/GANDRAS VISCI-G catheter deployment procedure should be performed by physicians thoroughly trained in percutaneous, intravascular techniques and procedures.

Do not use the GANDRAS/GANDRAS VISCI-G catheter if the packaging has been damaged.

Inspect the catheter prior to use for any kinks. Do not use a damaged catheter because vessel damage and/or inability to advance or withdraw the catheter may occur.

The catheter lumen must be flushed with sterile saline prior to use.

Precautions to prevent or reduce clotting should be taken when any catheter is used in the vascular system.

Exercise care while handling of the catheter during the procedure to reduce the possibility of accidental breakage, bending or kinking.

When the catheter is in the body, it should be manipulated only under fluoroscopy. Do not attempt to move the catheter without observing the resultant tip response.

If using the GANDRAS/GANDRAS VISCI-G catheter to deliver embolic materials, review any minimum catheter internal diameter recommendations specified by the maker of the therapeutic agent in order to minimize the potential for catheter obstruction. The minimum inner diameter of the GANDRAS/GANDRAS VISCI-G catheter is 0.037" (0.94mm).

CLINICAL PROCEDURE

The following instructions provide technical direction but do not obviate the necessity of formal training in the use of angiographic catheters. The techniques and procedures described do not represent ALL medically acceptable protocols, nor are they intended as a substitute for the physician's experience and judgment in treating any specific patient. All available data, including the patient's signs and symptoms and other diagnostic test results, should be considered before determining a specific treatment plan.

Package contains:

GANDRAS/GANDRAS VISCI-G catheter with attached tip straightener

Other materials required but not provided:

Sterile 0.035" (0.89mm) guidewire

Sterile syringe

Sterile heparinized saline

PREPARATIONS FOR USE

Carefully inspect the GANDRAS/GANDRAS VISCI-G catheter packaging and components for damage prior to use.

Utilizing sterile technique, remove the GANDRAS/GANDRAS VISCI-G catheter from its packaging and transfer to the sterile field.

Remove the catheter from packaging and inspect for any kinks. Remove the packaging stylet before use.

Thoroughly flush the catheter with heparinized saline.

Activate the hydrophilic coating by wetting the distal end of the catheter with sterile saline. For best results, keep the catheter surface wet at all times during use.

DEPLOYMENT PROCEDURE

Deploy the GANDRAS/GANDRAS VISCI-G catheter according to the following steps:

1. Load the inner lumen of a GANDRAS/GANDRAS VISCI-G catheter on a 0.035" (0.89mm) guidewire. The GANDRAS/GANDRAS VISCI-G catheter is supplied with a tip straightener to use at your discretion. To use the straightener, slide the straightener distally to open the catheter tip. With the catheter tip straight, advance the catheter over the wire. After initial advancement of the catheter over the wire, pull back and peel away the straightener. If not using the straightener, remove before use. The straightener must not be used as a vascular introducer sheath.
2. Insert the catheter through the already positioned introducer sheath of appropriate size.
3. Position the catheter in the proper location.

Warning: Never advance or withdraw an intravascular device against resistance until the cause of the resistance is determined by fluoroscopy. Movement of the catheter against resistance may result in separation of the catheter or guidewire tip, damage to the catheter or vessel perforation.

4. Remove the guidewire. Aspirate and flush the lumen and perform the angiographic procedure or deliver the appropriately prepared embolic materials.

Warning: Do not exceed the pressure rating listed on the label of the GANDRAS/GANDRAS VISCI-G catheter for any injection. Exceeding the pressure rating may result in separation of the catheter, embolization, vessel damage or delivery of the therapeutic agent to an unintended site.

PACKAGING & STORAGE

The GANDRAS/GANDRAS VISCI-G catheter has been sterilized with ethylene oxide.

STERILE EO

LIMITED WARRANTY

Vascular Solutions, Inc. warrants that the GANDRAS/GANDRAS VISCI-G catheter is free from defects in workmanship and materials prior to the stated expiration date. Liability under this warranty is limited to refund or replacement of any product, which has been found by Vascular Solutions, Inc. to be defective in workmanship or materials. Vascular Solutions, Inc. shall not be liable for any incidental, special or consequential damages arising from the use of the GANDRAS/GANDRAS VISCI-G catheter. Damage to the product through misuse, alteration, improper storage or improper handling shall void this limited warranty.

No employee, agent or distributor of Vascular Solutions, Inc. has any authority to alter or amend this limited warranty in any respect. Any purported alteration or amendment shall not be enforceable against Vascular Solutions, Inc.

THIS WARRANTY IS EXPRESSLY IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE OR ANY OTHER OBLIGATION OF VASCULAR SOLUTIONS, INC.

PATENTS AND TRADEMARKS

International and U.S. patents pending.

Gandras[®] is a registered trademark of Vascular Solutions, Inc.

See the International Symbols Glossary on page 34.

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