Minnie® Support Catheter

Minnie® Support Catheter

(Models 5700-5709) Instructions For Use

USA CAUTION

Federal (USA) law restricts this device to sale by or on the order of a physician.

CAUTION

The Minnie support catheters should be used by physicians with adequate training in the use of the device. The device is provided sterile for single use only. Reuse of single-use device creates a potential risk of patient or user infections. Contamination of the device may lead to illness or serious patient injury.

DEVICE DESCRIPTION

The Minnie support catheters are single lumen catheters designed for use in the arterial vasculature. The catheters provide support for guidewires during interventional procedures and allow for the exchange of one distally located guidewire for another one while maintaining access to distal vasculature. The Minnie support catheters have a hydrophilic coating on the distal portion of the catheter and are compatible with various guidewire diameters and lengths (see Table 1).

The Minnie catheters have three (3) evenly spaced radiopaque marker bands, the distal most of which is located 2.5mm from the tip of the catheter. Under fluoroscopy the radiopaque marker bands identify the distal tip of the catheter and can be used to estimate distances. Models 5700 and 5701 have white positioning marks located at 95cm (single mark) and 105cm (double marks) from the distal tip, respectively. The proximal end of the catheter incorporates a strain relief and a luer-lock guidewire entry port for flushing.

The Minnie support catheters are single lumen catheters designed to be used with guidewires with the following dimensions:

Table 1 – Specifications

Model Number	5709	5700	5701	5702
Max. guidewire diameter, in / mm	0.014" 0.36 mm	0.014" 0.36mm	0.014" 0.36mm	0.018" 0.46mm
Working length, cm "A"	90	135	150	90
Minimum guidewire length, cm	150	180	180	150
Marker band spacing, mm "B"	15	15	15	15
Marker band from distal tip, mm "C"	2.5	2.5	2.5	2.5
Proximal shaft diameter, in / mm	0.039" 0.99mm	0.039" 0.99mm	0.039" 0.99mm	0.044" 1.12mm
Distal shaft diameter, in / mm	0.026" 0.66mm	0.026" 0.66mm	0.026" 0.66mm	0.030" 0.76mm
Tip outside diameter, in / mm	0.020" 0.51mm	0.020" 0.51mm	0.020" 0.51mm	0.023" 0.58mm
Minimum Guide Catheter, F	5	5	5	6
Minimum Introducer Sheath, F	5	5	5	5


Model Number	5703	5704	5705
Max. guidewire diameter, in / mm	0.018" 0.46mm	0.018" 0.46mm	0.035" 0.89mm
Working length, cm "A"	135	150	65
Minimum guidewire length, cm	180	180	150
Marker band spacing, mm "B"	15	15	50
Marker band from distal tip, mm "C"	2.5	2.5	2.5
Proximal shaft diameter, in / mm	0.044" 1.12mm	0.044" 1.12mm	0.063" 1.60mm
Distal shaft diameter, in / mm	0.030" 0.76mm	0.030" 0.76mm	0.050" 1.27mm
Tip outside diameter, in / mm	0.023" 0.58mm	0.023" 0.58mm	0.041" 1.04mm
Minimum Guide Catheter, F	6	6	6
Minimum Introducer Sheath, F	5	5	5

Model Number	5706	5707	5708
Max. guidewire diameter, in / mm	0.035" 0.89mm	0.035" 0.89mm	0.035" 0.89mm
Working length, cm "A"	90	135	150
Minimum guidewire length, cm	150	180	180
Marker band spacing, mm "B"	50	50	50
Marker band from distal tip, mm "C"	2.5	2.5	2.5
Proximal shaft diameter, in / mm	0.063" 1.60mm	0.063" 1.60mm	0.063" 1.60mm
Distal shaft diameter, in / mm	0.050" 1.27mm	0.050" 1.27mm	0.050" 1.27mm
Tip outside diameter, in / mm	0.041" 1.04mm	0.041" 1.04mm	0.041" 1.04mm
Minimum Guide Catheter, F	6	6	6
Minimum Introducer Sheath, F	5	5	5

Table 2 – Minnie Catheter Flow Rates (ml/second) at 150psi and 300psi Injection Pressures for Saline and Contrast Solutions

Model	Size in/mm	Length cm	Sterile Saline		Contrast*	
			150 psi	300 psi	150 psi	300 psi
5709	0.014" 0.36mm	90	0.9	1.3	0.3	0.6
5700	0.014" 0.36mm	135	0.9	1.3	0.3	0.6
5701	0.014" 0.36mm	150	0.8	1.3	0.3	0.6
5702	0.018" 0.46mm	90	1.7	2.5	0.8	1.7
5703	0.018" 0.46mm	135	1.5	2.2	0.6	1.4
5704	0.018" 0.46mm	150	1.6	2.1	0.6	1.3
5705	0.035" 0.89mm	65	7.1	11.3	5.4	8.6
5706	0.035" 0.89mm	90	6.4	10.2	4.7	7.5
5707	0.035" 0.89mm	135	5.4	8.6	4.0	6.8
5708	0.035" 0.89mm	150	5.1	8.2	3.6	6.7

* Flow rates achieved with 4.3 cP contrast

INDICATIONS

The Minnie support catheters are intended to be used in conjunction with steerable guidewires in order to access discreet regions of the arterial and/or coronary vasculature. They may be used to facilitate placement and exchange of guidewires and other interventional devices. The Minnie support catheters also may be used to subselectively infuse/deliver therapeutic agents.

CONTRAINDICATIONS

The Minnie support catheters are contraindicated for use in synthetic vascular grafts due to the potential for graft perforation.

WARNINGS

Do not exceed maximum injection pressure of 300psi.

If an appropriate guidewire cannot be passed through the catheter, do not attempt to resolve the blockage by flushing the catheter *in vivo*. Catheter rupture and arterial injury could result. Identify and resolve the cause of the blockage or replace the catheter with a new one.

Never advance or withdraw an intravascular device against resistance until the cause of the resistance is determined by fluoroscopy. Movement of the catheter or guidewire against resistance may result in separation of the catheter or guidewire tip, damage to the catheter, or vessel perforation.

English/Instructions for Use	1
Český/Návod k použití	2
Dansk/Bruksanvisning	4
Nederlands/Gebruiksaanwijzing.....	5
Suomi/käyttöohjeet	6
Français/Mode d'emploi	8
Deutsch/Gebrauchsanweisung	9
Ελληνικά/Οδηγίες χρήσης	11
Italiano/Istruzioni per l'uso.....	12
Norsk/Bruksanvisning	14
Polski/Instrukcja użytkowania	15
Português/Instruções de Utilização.....	17
Español/Instrucciones de uso	18
Svenska/Bruksanvisning	19
Türkçe/Kullanım Talimatları	21

Vascular Solutions, Inc.
6464 Sycamore Court
Minneapolis, MN 55369 USA
(888) 240-6001 USA
(763) 656-4300
(763) 656-4250
www.vasc.com

Vascular Solutions Zerusa Limited
208 Business Innovation Centre
NUI Galway
Newcastle Road
Galway Ireland
+353 91 861611
+353 91 861612

COMPLICATIONS

As with all catheterization procedures, complications may occur when using the Minnie support catheters. These may include:

- local or systemic infection
- intimal disruption
- arterial dissection
- perforation and vessel rupture
- arterial thrombosis
- distal embolization of blood clots and plaque
- myocardial infarction
- arterial spasm
- catheter fracture with tip separation and distal embolization

PRECAUTIONS

The Minnie support catheter deployment procedure should be performed by physicians thoroughly trained in percutaneous, intravascular techniques and procedures.

Do not use the Minnie support catheter if the packaging has been damaged.

Inspect the catheter prior to use for any bends or kinks. Do not use a damaged catheter because vessel damage and/or inability to advance or withdraw the catheter may occur.

Care should be taken not to crush the catheter. Excessive tightening of a hemostatic valve onto the catheter shaft may result in damage to the wire lumen and difficulty while exchanging wires.

The entire catheter lumen should be thoroughly flushed with heparinized saline prior to use to prevent clot formation.

Precautions to prevent or reduce clotting should be taken when any catheter is used in the vascular system. Use of systemic heparinization and heparinized sterile solution should be considered.

Exercise care in handling the catheter during a procedure to reduce the possibility of accidental breakage, bending or kinking.

When the catheter is in the body, it should be manipulated only under fluoroscopy. Do not attempt to move the catheter without observing the resultant tip response.

CLINICAL PROCEDURE

The following instructions provide technical direction but do not obviate the necessity of formal training in the use of the Minnie catheters. The techniques and procedures described do not represent ALL medically acceptable protocols, nor are they intended as a substitute for the physician's experience and judgment in treating any specific patient.

Each Minnie support catheter includes the following components:

- Single-use disposable catheter
- Dispenser coil with flushing luer

Other materials required but not provided are:

- Catheter or sheath with an I.D. large enough to accommodate the proximal shaft O.D. of the specific model of Minnie support catheter in use (see Table 1)
- Guidewire with the appropriate O.D. and length for the model of Minnie catheter in use (see Table 1)
- 10ml syringe (for flushing the dispenser coil and the catheter lumen)
- Sterile heparinized saline (for system flushing)

PREPARATIONS FOR USE

1. Carefully inspect the Minnie catheter packaging and components for damage prior to use. Utilizing sterile technique, remove the Minnie catheter dispenser coil from its packaging and transfer it to the sterile field.
2. Attach a 10ml syringe filled with sterile, heparinized saline to the luer-lock guidewire entry port of the Minnie catheter and thoroughly flush the catheter.
3. Attach a 10ml syringe to the flushing luer on the dispenser coil and completely flush the coil to activate the catheter's hydrophilic coating.
4. Remove the Minnie catheter from the dispenser coil and inspect it for any bends or kinks.

DEPLOYMENT PROCEDURE

The following Minnie catheter deployment steps assume a standard catheterization protocol for exchange of a guidewire.

As with any interventional procedure, proper

anticoagulation and anti-platelet therapy should be administered prior to beginning.

Note: Familiarity with traditional guidewire exchange techniques is required for successful deployment of the Minnie catheter and guidewire.

When not in use during this procedure, wipe the Minnie catheter with a sterile gauze pad saturated with heparinized saline, flush the lumen well, and store in the dispenser coil in a saline bath.

1. Backload the distal tip of the Minnie catheter onto the proximal end of the appropriate length guidewire that is already in place in the distal vasculature.
 2. Carefully advance the catheter until the marker band at the distal tip is visible in the desired distal vascular space.
- WARNING: Never advance or withdraw an intravascular device against resistance until the cause of the resistance is determined by fluoroscopy. Movement of the catheter or guidewire against resistance may result in separation of the catheter or guidewire tip, damage to the catheter, or vessel perforation.**
3. After fluoroscopically confirming that the catheter tip is positioned in the desired vascular space, slowly remove the guidewire.
 4. Insert the desired length replacement guidewire into the luer-lock guidewire entry port. Advance the guidewire until it exits the catheter tip.
 5. Under fluoroscopy, ensure that the guidewire is sufficiently advanced into the distal vasculature. Fix the guidewire using standard guidewire exchange techniques and carefully withdraw the catheter.

PACKAGING & STORAGE

The Minnie support catheter has been sterilized with ethylene oxide.

STERILE EO

Handle with care.

LIMITED WARRANTY

Vascular Solutions, Inc. warrants that the Minnie support catheter is free from defects in workmanship and materials prior to the stated expiration date. Liability under this warranty is limited to refund or replacement of any product which has been found by Vascular Solutions, Inc. to be defective in workmanship or materials. Vascular Solutions, Inc. shall not be liable for any incidental, special or consequential damages arising from the use of the MINNIE support catheter. Damage to the product through misuse, alteration, improper storage or improper handling shall void this limited warranty.

No employee, agent or distributor of Vascular Solutions, Inc. has any authority to alter or amend this limited warranty in any respect. Any purported alteration or amendment shall not be enforceable against Vascular Solutions, Inc.

THIS WARRANTY IS EXPRESSLY IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE OR ANY OTHER OBLIGATION OF VASCULAR SOLUTIONS, INC.

PATENTS AND TRADEMARKS

International and U.S. patents pending.

Minnie® is a registered trademark of Vascular Solutions, Inc.

See the International Symbols Glossary on page 23.

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Nigdy nie przesuwaj i nie wycofuj urządzenia wewnątrznaczyniowego po napotkaniu oporu, dopóki przyczyna oporu nie zostanie ustalona fluoroskopowo. Przesuwanie cewnika lub przewodnika w przypadku wystąpienia oporu może spowodować oddzielenie się końcówki cewnika lub © 2011 Vascular Solutions, Inc.
Minnie® Vascular Solutions, Inc.'in tescilli bir ticari markasıdır.