English/Instructions for Use	1
Español/Instrucciones de uso.....	3

Vascular Solutions, Inc.
6464 Sycamore Court
Minneapolis, MN 55369
USA
(888) 240-6001 USA
(763) 656-4300
(763) 656-4250
www.vascularsolutions.com

D-Stat® Dry Wrap Hemostatic Bandage Instructions For Use

USA CAUTION

Federal (USA) law restricts this device to sale by or on the order of a physician.

CAUTION

The D-Stat Dry Wrap procedure should be performed by physicians or physician directed allied health care professionals with adequate training in the use of the device.

DEVICE DESCRIPTION

Each D-Stat Dry Wrap hemostatic bandage consists of a lyophilized pad consisting of thrombin, sodium carboxymethylcellulose and calcium chloride in a nonwoven gauze.

Thrombin is a protein substance produced through a conversion reaction in which prothrombin of bovine origin is activated by tissue thromboplastin of bovine origin in the presence of calcium chloride. Thrombin contains no preservative and has been chromatographically purified. Thrombin requires no intermediate physiological agent for its reaction. It converts fibrinogen directly to fibrin. This product contains not less than 200 units of bovine-derived thrombin.

Sodium carboxymethylcellulose, also known as cellulose gum or CMC, serves as the matrix for the lyophilized pad and as a suspension agent for the thrombin. Calcium chloride is added to assist in the clotting cascade.

Hemostasis is achieved by the physiological coagulation-inducing properties of the lyophilized pad combined with manual compression.

The D-Stat Dry Wrap has been sterilized with irradiation.

INDICATIONS

The D-Stat Dry Wrap hemostatic bandage is applied topically as an adjunct to manual compression and is indicated for the control of surface bleeding from vascular access sites and percutaneous catheters or tubes.

CONTRAINDICATIONS

The D-Stat Dry Wrap is contraindicated in persons with known sensitivity to bovine-derived materials.

WARNINGS

The use of topical bovine thrombin preparations has occasionally been associated with abnormalities in hemostasis ranging from asymptomatic alterations in laboratory determinations, such as prothrombin time (PT) and partial thromboplastin time (PTT), to severe bleeding or thrombosis which rarely have been fatal. These hemostatic effects appear to be related to the formation of antibodies against bovine thrombin and/or factor V which in some cases may cross react with human factor V, potentially resulting in factor V deficiency. Repeated clinical applications of topical bovine thrombin increase the likelihood that antibodies against thrombin and/or factor V may be formed. Consultation with an expert in coagulation disorders is recommended if a patient exhibits abnormal coagulation laboratory values, abnormal bleeding, or abnormal thrombosis following the use of topical thrombin. Any interventions should consider the immunologic basis of this condition. Patients with antibodies to bovine thrombin preparations should not be re-exposed to these products.

The D-Stat Dry Wrap is supplied sterile for single use only. Do not re-sterilize.

STERILE | R

Do not use the D-Stat Dry Wrap as a replacement for absorbable hemostats. This product contains non-absorbable materials and is not intended to be left in the body.

Do not place the D-Stat Dry Wrap into blood vessels. Extensive intravascular clotting and even death may result.

PRECAUTIONS

Do not use the D-Stat Dry Wrap if the packaging has been damaged.

The safety and effectiveness of the D-Stat Dry Wrap have not been established in children and pregnant women.

Do not touch the D-Stat Dry Wrap lyophilized pad with wet gloves or expose the pad to liquid. Absorption of the liquid and destruction of the pad will result.

The D-Stat Dry Wrap should not be used in the presence of infection. It should be used with caution in contaminated areas of the body.

Do not leave the D-Stat Dry Wrap adhesive bandage applied for more than 24 hours. Skin irritation may result.

When removing the D-Stat Dry Wrap lyophilized pad and the adhesive bandage, do not disrupt the clot by physical manipulation. If the pad adheres to the puncture site, irrigate the pad with non-heparinized saline and carefully remove it.

ADVERSE EVENTS

A recognized rare potential reaction associated with the use of bovine-derived thrombin is the development of inhibitory antibodies, which interferes with hemostasis.

CLINICAL PROCEDURE

The following instructions provide technical direction but do not obviate the necessity of formal training in the use of the D-Stat Dry Wrap hemostatic bandage. The techniques and procedures described do not represent ALL medically acceptable protocols, nor are they intended as a substitute for the experience and judgment of the physician or trained clinical professional in treating any specific patient.

Carefully inspect the D-Stat Dry Wrap hemostatic bandage packaging for damage prior to use.

APPLICATION PROCEDURE

- Using sterile technique, open the foil pouch and remove the D-Stat Dry Wrap pad. Transfer the D-Stat Dry Wrap into the sterile field.
- Place the D-Stat Dry Wrap directly over the source of the bleeding and apply adjunct manual pressure. The pre-formed slit in the D-Stat Dry Wrap allows for placement around indwelling lines, if desired.
- Once hemostasis has been achieved, leave the D-Stat Dry Wrap in place and apply adhesive bandage.
- If desired, the D-Stat Dry Wrap may be left in place for up to 24 hours.

Precaution: Do not leave the D-Stat Dry Wrap or adhesive bandage applied for more than 24 hours. Skin irritation may result.

Precaution: Upon removal, do not disrupt the clot by physical manipulation. If the D-Stat Dry Wrap or adhesive bandage adheres to the placement site, gently irrigate the pad and bandage with non-heparinized saline and carefully remove it.

3M™ Tegaderm™ Transparent Film Dressing Frame Style

1626, 1627, 1628, 1629, 1630, 1634, 1622W, 1623W, 1624W, 1626W, 9505W, 9506W

Description: Tegaderm™ Film consists of a thin film backing with a non-latex, hypoallergenic adhesive. Tegaderm™ Film with Border is notched and reinforced with soft cloth tape to provide a better seal around catheters and other devices. The dressing is breathable, allowing good oxygen and moisture vapor exchange. It is waterproof and impermeable to liquids, bacteria, and viruses. *An intact dressing protects the site from outside contamination.


*Laboratory testing has proven Tegaderm™ and Tegaderm™ HP Film provide a barrier against HIV-1 and HBV while the dressings remain intact without leakage.

Indications: Tegaderm™ Film can be used to cover and protect catheter sites and wounds, to maintain a moist environment for wound healing or to facilitate autolytic debridement, as a secondary dressing, as a protective cover over at-risk skin, to secure devices to the skin, to cover first and second degree burns, and as a protective eye covering. Do not use the dressing as a replacement for sutures and other primary wound closure methods.

Precautions:

1. Stop any bleeding at the site before applying the dressing.
2. Do not stretch the dressing during application as tension can cause skin trauma.
3. Make sure the skin is clean, free of soap residue and lotion and allowed to dry thoroughly before applying the dressing to prevent skin irritation and to ensure good adhesion.
4. The dressing may be used on an infected site, only when under the care of a health care professional.
5. Antimicrobial ointments containing polyethylene glycols may compromise the strength of the Tegaderm™ HP Transparent Film Dressings.
6. Tegaderm™ Transparent Dressings should not be re-sterilized by gamma, E-beam or steam methods.

Instructions for Use: Refer to figures.

If you have any questions or comments, contact the 3M Health Care Customer Help Line at **1-800-228-3957 OR GO TO WWW.3M.COM.**

LIMITED WARRANTY

Vascular Solutions, Inc. warrants that the D-Stat Dry Wrap hemostatic bandage is free from defects in workmanship and materials prior to the stated expiration date. Liability under this warranty is limited to refund or replacement of any product which has been found by Vascular Solutions, Inc. to be defective in workmanship or materials. Vascular Solutions, Inc. shall not be liable for any incidental, special, or consequential damages arising from the use of D-Stat Dry Wrap hemostatic bandage. Damage to the product through misuse, alteration, improper storage, or improper handling shall void this limited warranty.

No employee, agent, or distributor of Vascular Solutions, Inc. has any authority to alter or amend this limited warranty in any respect. Any purported alteration or amendment shall not be enforceable against Vascular Solutions, Inc.

THIS WARRANTY IS EXPRESSLY IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING ANY THIS WARRANTY IS EXPRESSLY IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE OR ANY OTHER OBLIGATION OF VASCULAR SOLUTIONS, INC.

PATENTS AND TRADEMARKS

U.S. patents pending.

D-Stat® is a registered trademark of Vascular Solutions, Inc.

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D-Stat[®] es una marca registrada de Vascular Solutions, Inc.