

PRODUCT LAUNCH ANNOUNCEMENT



Tuesday, March 3, 2009

Contacts: Bob Oas, Senior Marketing Communications Specialist, boas@vascularsolutions.com
Joan Will, Director of Marketing Communications, jwill@vascularsolutions.com
Vascular Solutions, Inc. 763.656.4300

Vascular Solutions Launches Flamingo™ Inflation Device and Interventional Packs

MINNEAPOLIS, Minnesota -- Vascular Solutions, Inc. recently announced its partnership with Sedat and the launch of the Flamingo inflation device. The Flamingo is designed for use during cardiovascular procedures to create, maintain and monitor pressures up to 30atm in a compact 20cc syringe. The Flamingo offers progressive inflation resistance while providing constant pressure feedback during inflation. It also features an ergonomic design with a 90° rotating gauge for customized handling, and it is constructed of clear materials to ensure visualization of air bubbles.



The Flamingo inflation device comes with an optional 3-way stopcock and Alligator™ torque device. Vascular Solutions also offers convenient interventional packs complete with a Flamingo inflation device, a Guardian® hemostasis valve and other optional accessories. Developed and manufactured by Sedat, the Flamingo is currently available in the United States through Vascular Solutions' sales force.

Vascular Solutions, Inc. is a leading medical device company that delivers proprietary clinical solutions for diagnosing and treating vascular conditions. The company's rapidly growing product line consists of innovative devices across established and emerging areas of coronary and peripheral vascular medicine. For more information, visit www.vascularsolutions.com.

Flamingo™ and Alligator™ are trademarks of Sedat. Guardian® is a registered trademark of Zerusa Ltd.