

PRODUCT LAUNCH ANNOUNCEMENT



Friday, March 6, 2009

Contacts: Bob Oas, Senior Marketing Communications Specialist, boas@vascularsolutions.com
Joan Will, Director of Marketing Communications, jwill@vascularsolutions.com
Vascular Solutions, Inc. 763.656.4300

Vascular Solutions Launches Unique GrebSet™ Micro-Introducer Kit

MINNEAPOLIS, Minnesota -- Vascular Solutions, Inc. recently announced the launch of its GrebSet, a unique micro-introducer kit with an angled tip that facilitates steerability.



The GrebSet is a 30cm coaxial micro-introducer kit that allows the operator to direct an .035" or .038" guidewire into a vessel branch without having to exchange from a conventional access kit to an angled catheter after gaining access. Additional features include hydrophilic coating for smooth insertion, a braided sheath for excellent torque, and a radiopaque marker band embedded in the distal tip for enhanced visualization. Currently available in the U.S., the GrebSet comes packaged with an .018" nitinol guidewire and either a 21G 15cm Trocar needle or 7cm single-wall echogenic needle.

Vascular Solutions, Inc. is a leading medical device company that delivers proprietary clinical solutions for diagnosing and treating vascular conditions. The company's rapidly growing product line consists of innovative devices across established and emerging areas of coronary and peripheral vascular medicine. For more information, visit www.vascularsolutions.com.