

PRODUCT LAUNCH ANNOUNCEMENT



Thursday, July 24, 2008

Contacts: Bob Oas, Senior Marketing Communications Specialist, boas@vascularsolutions.com
Joan Will, Director of Marketing Communications, jwill@vascularsolutions.com
Vascular Solutions, Inc. 763.656.4300

Vascular Solutions Launches Jiffy™ 0.035" Specialty Guidewire

MINNEAPOLIS, Minnesota — Vascular Solutions, Inc. recently launched the Jiffy Wire, an 80cm 0.035" specialty guidewire.

The Jiffy Wire offers the benefits of both a straight floppy tip and a 3mm J-tip in one double-ended guidewire. The 80cm working length and stiff shaft of the Jiffy Wire provides control and eases wire manipulation for access procedures, while the 10cm graduated transition tips facilitate atraumatic access and device placement. The coils of the Jiffy Wire are pre-coated with a smooth PTFE coating. The Jiffy Wire is currently available in the United States.



Vascular Solutions, Inc. is a leading medical device company that delivers proprietary clinical solutions for diagnosing and treating vascular conditions. The company's rapidly growing product line consists of innovative devices across established and emerging areas of coronary and peripheral vascular medicine. For more information, visit www.vascularsolutions.com.