


Gandras[®]

Catheters

New advances in catheterization of uterine arteries and visceral vascular territories



- Ultra-tapered distal segment with hydrophilic coating for the benefits of a microcatheter without the added cost
- Available in 2cm, 4cm and 6cm distal tip lengths for customized arterial placement—NOW with marker band at primary curve (models 5580, 5581, 5582)
- Flexible shaft facilitates ipsilateral and contralateral access
- Radiopaque distal segment for enhanced visualization
- NEW VISCI-G (model 5585) combines shorter overall length and shorter distal tip to facilitate precise and stable arterial placement in the visceral vascular territories

Gandras[®] Catheter


Gandras[®] Visci-G Catheter

The Gandras/Gandras Visci-G catheter are hydrophilically coated catheters for use in delivering embolic materials and radiopaque media to selected sites in the vascular system. The choice of a specific catheter for an angiographic procedure should be based on the physician's judgment and experience.



Model	Length	Tip Length	Platinum-Iridium Marker Band	Proximal O.D.	Distal O.D.	Recommended Guidewire	Pressure Limit	Box Quantity*
5580	82cm	2cm	Yes	5F	4F	0.035"	600psi	5
5581	84cm	4cm	Yes	5F	4F	0.035"	600psi	5
5582	86cm	6cm	Yes	5F	4F	0.035"	600psi	5
5585	75cm	1cm	No	5F	5F	0.035"	600psi	5

* Devices packaged in quantities of 5 units per box. Individual units are not sold separately.



The Gandras/Gandras Visci-G catheter is designed to be used for delivering embolic materials and radiopaque media to selected sites in the vascular system. Diagnostic, embolic, or therapeutic agents are to be used in accordance with specifications outlined by the manufacturer. Please see the Instructions for Use for a complete listing of indications, contraindications, warnings and precautions.

CAUTION: Federal law (U.S.A.) restricts these devices to sale by or on the order of a physician.

Gandras is a registered trademark of Vascular Solutions, Inc. U.S. and International patents pending.
©2010 Vascular Solutions, Inc. All rights reserved.
ML1879 Rev. F 07/10



Vascular Solutions, Inc.
6464 Sycamore Court
Minneapolis, Minnesota 55369 USA

Customer Service:
United States: 888.240.6001
International: (001) 763.656.4298
customerservice@vascularsolutions.com



www.vascularsolutions.com