

PRONTO[®]

Extraction Catheters

Think
PRONTO[®]
Extraction Catheters
FIRST

The choice of aspiration catheters.

	PRONTO V₃ Extraction Catheter	PRONTO LP Extraction Catheter	PRONTO .035" Extraction Catheter	PRONTO-SHORT Extraction Catheter
Vessel Size	≥2.0 mm	≥1.5mm	≥4.0mm	≥2.5mm
Soft Thrombus Removal From:*	Coronary and Peripheral Arteries	Coronary and Peripheral Arteries	Large Vessels, Veins and Arteries	Synthetic Dialysis Grafts and AV Fistulas
Guide/Sheath Compatibility	6F Guide (0.070" I.D.)	6F Guide (0.066" I.D.)	10F Sheath	6F Sheath
Crossing Profile	0.065"	0.053"	0.128"	0.078"
Guidewire Compatibility	0.014"	0.014"	0.035"	0.018"
Guidewire Lumen	Rapid Exchange	Rapid Exchange	OTW	OTW
Working Length	140cm	140cm	115cm	40cm
Distal Coating	Hydrophilic	Hydrophilic	Siliconized	Bare
Model Number	5003	5010	5040	5030

*The Pronto Family includes catheter models indicated for the removal of fresh, soft emboli and thrombi from coronary arteries, veins, bypass grafts, AV access conduits, and/or peripheral arteries. Please see each Pronto product Instructions for Use for a complete listing of the indications, contraindications, warnings and precautions.

U.S.A. CAUTION: Federal Law (U.S.A.) restricts this device to sale by or on the order of a physician.
CAUTION: The Pronto Extraction Catheter procedures should be performed by physicians with adequate training in the use of the device. Pronto is a registered trademark of Vascular Solutions, Inc.
©2008 Vascular Solutions, Inc. All rights reserved. ML1917 Rev. B 11/09

Vascular Solutions, Inc.
6464 Sycamore Court
Minneapolis, Minnesota 55369 USA
USA Customer Service: 888.240.6001
www.vascularsolutions.com



vascular
SOLUTIONS