


Jiffy™ Guidewire



The Jiffy Wire is an 0.035" guidewire that offers the benefits of both a straight floppy tip and a 3mm J-tip all in one guidewire.

- Double-ended tip design features advantage of two guidewires in one
- 10cm graduated tips facilitate atraumatic access and device placement
- 80-150cm working length provides control and eases wire manipulation for access procedures
- Peel-away J-tip straightener included for convenient delivery

Jiffy™ Guidewire

The Jiffy Wire is intended to facilitate the placement of devices during diagnostic and interventional procedures.

Model	Wire Diameter	Wire Length	Tip	Box Quantity*
7380	0.035"	80cm	Double-ended, Straight Floppy & 3mm J	5
7381	0.035"	120cm	Double-ended, Straight Floppy & 3mm J	5
7382	0.035"	150cm	Double-ended, Straight Floppy & 3mm J	5

*Devices packaged in quantities of 5 units per box. Individual units are not sold separately.




VSI Torque Device accessory

The VSI Torque Device is designed to be used with guidewires ranging in size from 0.014" to 0.038".

Model	Description	Box Quantity*
8180	VSI Torque Device	10

*Devices packaged in quantities of 10 units per box. Individual units are not sold separately.

The VSI Torque Device is intended for use in manipulating a guidewire that is in the vasculature.



vascular
SOLUTIONS

Manufactured for:
Vascular Solutions, Inc.
 6464 Sycamore Court
 Minneapolis, Minnesota 55369
 Customer Service: 888.240.6001
 www.vasc.com

Please see the Instructions for Use for a complete listing of the indications, contraindications, warnings and precautions.

CAUTION: Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.

Jiffy is a trademark of Vascular Solutions, Inc.

©2011 Vascular Solutions, Inc. All rights reserved. ML1928 Rev. D 12/11