

Amplatz **SST**[™]

Guidewire

*Super stiff core 0.035" and 0.038"
guidewires for extra support*

- Inner-core construction provides added strength for increased support
- Radiopaque tip for enhanced visibility
- PTFE coating for enhanced device tracking
- Available in straight and modified J configurations

Amplatz **SST**[™]

Guidewire

The Amplatz SST steerable guidewire facilitates placement of a catheter during diagnostic angiography and interventional procedures. The Amplatz SST can be accurately controlled to facilitate navigation through tortuous vessels and/or adjoining side branches.

Amplatz SST, Straight

Model	Wire Diameter	Wire Length	Tip Shape	Box Quantity*
7479	0.035"	80cm	Straight	5
7480	0.035"	150cm	Straight	5
7481	0.035"	180cm	Straight	5
7482	0.035"	260cm	Straight	5
7483	0.038"	80cm	Straight	5
7484	0.038"	150cm	Straight	5
7485	0.038"	180cm	Straight	5
7486	0.038"	260cm	Straight	5

Amplatz SST, Modified J

Model	Wire Diameter	Wire Length	Tip Shape	Box Quantity*
7487	0.035"	150cm	Modified J	5
7488	0.035"	180cm	Modified J	5
7489	0.035"	260cm	Modified J	5
7490	0.038"	150cm	Modified J	5
7491	0.038"	180cm	Modified J	5

*Packaged in quantities of 5 units per box. Individual units are not sold separately.

VSI Torque Device accessory

The VSI Torque Device is designed to be used with guidewires ranging in size from 0.014" to 0.038".

Model	Description	Box Quantity*
8180	VSI Torque Device	10

*Devices packaged in quantities of 10 units per box. Individual units are not sold separately.

The VSI Torque Device is intended for use in manipulating a guidewire that is in the vasculature.




Please see the Instructions for Use for a complete listing of the indications, contraindications, warnings and precautions.

CAUTION: Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.

Amplatz SST is a trademark of Vascular Solutions, Inc.

©2010 Vascular Solutions, Inc. All rights reserved. ML2039 Rev. B 07/10



vascular
SOLUTIONS

Manufactured for:
Vascular Solutions, Inc.

6464 Sycamore Court
Minneapolis, Minnesota 55369
Customer Service: 888.240.6001
www.vascularsolutions.com