

vsi Tru-Torque™

Guidewire

A complete line of highly torqueable 0.035" guidewires

- Excellent torque response for steerability
- Radiopaque tip for enhanced visibility
- Packaged with a torque device
- Available in standard, floppy and modified J configurations

VSI Tru-Torque™

Guidewire

The Tru-Torque steerable guidewire facilitates placement of a catheter during diagnostic angiography and interventional procedures. The Tru-Torque can be accurately controlled to facilitate navigation through tortuous vessels and/or adjoining side branches. A torque device is included with each Tru-Torque guidewire.

VSI Tru-Torque, Standard—highly supportive

Model	Wire Diameter	Wire Length	Tip Shape	Box Quantity*
7470	0.035"	150cm	Straight	5
7471	0.035"	300cm	Straight	5

VSI Tru-Torque, Floppy—highly flexible

Model	Wire Diameter	Wire Length	Tip Shape	Box Quantity*
7472	0.035"	150cm	Straight	5
7473	0.035"	180cm	Straight	5
7474	0.035"	260cm	Straight	5

VSI Tru-Torque, Modified J—highly versatile

Model	Wire Diameter	Wire Length	Tip Shape	Box Quantity*
7475	0.035"	150cm	Modified J	5
7476	0.035"	180cm	Modified J	5
7477	0.035"	260cm	Modified J	5
7478	0.035"	300cm	Modified J	5

VSI Tru-Torque, Modified J with Extra Flop**—highly versatile

Model	Wire Diameter	Wire Length	Tip Shape	Box Quantity*
7469	0.035"	150cm	Modified J	5

*Packaged in quantities of 5 units per box. Individual units are not sold separately.

**Extra Flop model contains an additional 2cm of flexibility on the distal floppy tip.

VSI Torque Device accessory

The VSI Torque Device is designed to be used with guidewires ranging in size from 0.014" to 0.038".

Model	Description	Box Quantity*
8180	VSI Torque Device	10

*Devices packaged in quantities of 10 units per box. Individual units are not sold separately.


The VSI Torque Device is intended for use in manipulating a guidewire that is in the vasculature.

Please see the Instructions for Use for a complete listing of the indications, contraindications, warnings and precautions.

CAUTION: Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.

VSI Tru-Torque is a trademark of Vascular Solutions, Inc.

©2010 Vascular Solutions, Inc. All rights reserved. ML2040 Rev. C 07/10



Vascular
SOLUTIONS

Manufactured for:
Vascular Solutions, Inc.
 6464 Sycamore Court
 Minneapolis, Minnesota 55369
 Customer Service: 888.240.6001
www.vascularsolutions.com