



Rad-Band™ *Topical Hemostat*

***Easy-to-use radial
compression band***



Hands-free pressure in a comfortable hemostasis band

- Easy release clasp for single operator adjustment of pressure
- Clear compression window allows for continuous site monitoring
- Adjustable foam pads for patient's comfort

Rad-Band™

Topical Hemostat

Model	Description	Box Quantity
3505	Rad-Band Topical Hemostat	10*

* Devices packaged in quantities of 10 units per box. Individual units are not sold separately.

The Rad-Band topical hemostat is applied topically and is indicated for the control of surface bleeding from vascular access sites and percutaneous catheters or tubes.



Please see the Instructions for Use for a complete listing of the indications, contraindications, warnings and precautions.
CAUTION: Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.

U.S. and International patents pending.
Rad-Band is a trademark of Vascular Solutions, Inc.

©2011 Vascular Solutions, Inc. All rights reserved.
ML2126 Rev. B 01/11



Vascular Solutions, Inc.
6464 Sycamore Court
Minneapolis, Minnesota 55369 USA
Customer Service: (888) 240-6001
www.vascularsolutions.com