



D-Stat[®] Radial

*Hemostatic Band
Topical Hemostat*

*Utilizes the power of
thrombin to facilitate
hemostasis*

More than just a radial compression band

- Contains thrombin in familiar D-Stat Dry pad
- Designed to prevent compression of the ulnar artery
- Low-profile band allows wrist full range-of-motion
- Easy to tighten and loosen without loss of pressure
- Adjustable foam pads for patient's comfort

D-Stat[®] Radial

Hemostatic Band

Topical Hemostat

| Model | Description | Box Quantity |
|-------|-----------------------------|--------------|
| 3501 | D-Stat Radial Hemostat Band | 10* |

* Devices packaged in quantities of 10 units per box. Individual units are not sold separately.

The D-Stat Radial hemostat band is indicated for the local management and control of bleeding from vascular access sites and percutaneous catheters and tubes.



The use of topical bovine thrombin preparations has occasionally been associated with abnormalities in hemostasis ranging from asymptomatic alterations in laboratory determinations, such as prothrombin time (PT) and partial thromboplastin time (PTT), to severe bleeding or thrombosis which rarely have been fatal. These hemostatic effects appear to be related to the formation of antibodies against bovine thrombin and/or factor V which in some cases may cross react with human factor V, potentially resulting in factor V deficiency. Repeated clinical applications of topical bovine thrombin increase the likelihood that antibodies against thrombin and/or factor V may be formed. Consultation with an expert in coagulation disorders is recommended if a patient exhibits abnormal coagulation laboratory values, abnormal bleeding, or abnormal thrombosis following the use of topical thrombin. Any interventions should consider the immunologic basis of this condition. Patients with antibodies to bovine thrombin preparations should not be re-exposed to these products.

The D-Stat Radial is contraindicated in persons with known sensitivity to bovine-derived materials. Please see the Instructions for Use for a complete listing of the indications, contraindications, warnings and precautions.

CAUTION: Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.

CAUTION: The D-Stat Radial procedure should be performed by physicians or physician-directed allied health care professionals with adequate training in the use of the device.

U.S. and International patents pending.
D-Stat is a registered trademark of Vascular Solutions, Inc.

©2009 Vascular Solutions, Inc. All rights reserved.
ML2173 Rev. A 12/09



Vascular Solutions, Inc.
6464 Sycamore Court
Minneapolis, Minnesota 55369 USA

Customer Service:
United States: 888.240.6001
International: (001) 763.656.4298
customerservice@vascularsolutions.com



www.vascularsolutions.com