

PRODUCT LAUNCH ANNOUNCEMENT




Monday, May 4, 2009

Contacts: Joan Will, Director of Marketing Communications, jwill@vascularsolutions.com
Vascular Solutions, Inc. 763.656.4300

Vascular Solutions Launches VSI Tru-Torque™ guidewires

MINNEAPOLIS, Minnesota -- Vascular Solutions, Inc. recently launched the VSI Tru-Torque line of guidewires, a complete line of highly torqueable 0.035" guidewires that facilitate catheter placement during diagnostic angiography and interventional procedures.



The VSI Tru-Torque guidewires feature an excellent torque response for steerability, PTFE coating for enhanced device tracking and radiopaque tip for enhanced visibility. Each VSI Tru-Torque guidewire comes packaged with a torque device. The VSI Tru-Torque guidewires can be accurately controlled to facilitate navigation through tortuous vessels and adjoining side branches. They are available in standard, floppy and modified J tip configurations, and are currently available in the United States.

Vascular Solutions, Inc. is a leading medical device company that delivers proprietary clinical solutions for diagnosing and treating vascular conditions. The company's rapidly growing product line consists of innovative devices across established and emerging areas of coronary and peripheral vascular medicine. For more information, visit www.vascularsolutions.com.