

PRODUCT LAUNCH ANNOUNCEMENT




Wednesday, January 20, 2010

Contacts: Joan Will, Director of Marketing Communications, jwill@vascularsolutions.com
Vascular Solutions, Inc. 763.656.4300

Vascular Solutions Expands Jiffy™ 0.035" Specialty Guidewire Line

MINNEAPOLIS, Minnesota -- Vascular Solutions, Inc. recently introduced two new Jiffy™ guidewires, adding both a 120cm and 150cm length to the current 80cm, 0.035" specialty guidewire line.

The Jiffy guidewire offers the benefits of both a straight floppy tip and a 3mm J-tip in one double-ended guidewire. The stiff shaft of the Jiffy guidewire provides control and eases wire manipulation for access procedures, while the 10cm graduated transition tip facilitates atraumatic access and device placement. A smooth PTFE pre-coating provides added lubricity. The Jiffy guidewire is currently available in the United States.



Vascular Solutions, Inc. is a leading medical device company that delivers proprietary clinical solutions for diagnosing and treating vascular conditions. The company's rapidly growing product line consists of innovative devices across established and emerging areas of coronary and peripheral vascular medicine. For more information, visit www.vascularsolutions.com.