

# PRODUCT LAUNCH ANNOUNCEMENT



Monday, February 1, 2010

Contacts: Joan Will, Director of Marketing Communications, [jwill@vascularsolutions.com](mailto:jwill@vascularsolutions.com)  
Vascular Solutions, Inc. 763.656.4300

## Vascular Solutions Launches Gator™ ClipSeal Plug

MINNEAPOLIS, Minnesota -- Vascular Solutions, Inc. announces the availability of the new Gator™ ClipSeal Plug designed to maintain a hemostatic seal around 0.035" or 0.038" guidewires within large hemostatic valve introducer sheaths to minimize blood loss. The Gator is intended for use with commonly used 12F to 24F introducer sheaths in procedures such as Endovascular Aneurysm Repair (EVAR) to treat abdominal aortic aneurysms (AAA) or thoracic aortic aneurysms (TAA).



The unique design of the Gator provides an easy and effective solution for minimizing blood loss by side-loading onto the in-place guidewire and securely fitting into the valve of the introducer sheath. The Gator is currently available for sale in the United States through Vascular Solutions' direct sales force.

Vascular Solutions, Inc. is a leading medical device company that delivers proprietary clinical solutions for diagnosing and treating vascular conditions. The company's rapidly growing product line consists of innovative devices across established and emerging areas of coronary and peripheral vascular medicine. For more information, visit [www.vascularsolutions.com](http://www.vascularsolutions.com).