

# PRODUCT LAUNCH ANNOUNCEMENT



Monday, April 12, 2010

Contacts: Joan Will, Director of Marketing Communications, [jwill@vascularsolutions.com](mailto:jwill@vascularsolutions.com)  
Vascular Solutions, Inc. 763.656.4300

## Vascular Solutions Launches New Gandras® Visci-G Catheter

MINNEAPOLIS, Minnesota -- Vascular Solutions, Inc. recently launched the Gandras® Visci-G catheter for use in delivering embolic materials and radiopaque media to selected sites in visceral vascular territories.

The Gandras Visci-G is a 5F, 75 cm catheter with a unique primary curve and a flexible 1 cm tip, making it useful for precise and stable placement in renal, mesenteric or abdominal catheterization procedures. The tight primary curve allows easy access of bilateral arteries through a single 5F arterial introducer, and standard microcatheters can then be delivered through the Visci-G catheter to access distal regions. The Gandras Visci-G catheter can be used in clinical cases such as GI hemorrhages, endograft leaks, and anytime operators are having difficulty seating a catheter in an artery that branches off the aorta. The new Gandras Visci-G is now available worldwide through Vascular Solutions' clinically trained U.S. sales force and international distributors.



Vascular Solutions, Inc. is a leading medical device company that delivers proprietary clinical solutions for diagnosing and treating vascular conditions. The company's rapidly growing product line consists of innovative devices across established and emerging areas of coronary and peripheral vascular medicine. For more information, visit [www.vascularsolutions.com](http://www.vascularsolutions.com).