

# PRODUCT LAUNCH ANNOUNCEMENT




Wednesday, March 19, 2008

Contacts: Bob Oas, Senior Marketing Communications Specialist, [boas@vascularsolutions.com](mailto:boas@vascularsolutions.com)  
Joan Will, Director of Marketing Communications, [jwill@vascularsolutions.com](mailto:jwill@vascularsolutions.com)  
Vascular Solutions, Inc. 763.656.4300

## Vascular Solutions Launches New Siliconized Vari-Lase® Introducer Sheaths for Use in Endovenous Laser Therapy

MINNEAPOLIS, Minnesota — Vascular Solutions, Inc. recently added to its Vari-Lase endovenous laser therapy offering with the introduction of new siliconized introducer sheaths. The new sheaths allow smooth insertion into and easy, even pull-back from the vein during the treatment of varicose veins with endovenous laser energy.



Now included in all Vari-Lase Standard and Completely Custom endovenous laser therapy procedure kits, the new siliconized sheaths provide an approximated 50% reduction in pull-back force compared with conventional untreated introducer sheaths. The siliconized Vari-Lase introducer sheaths are currently available worldwide.

Vascular Solutions, Inc. is a leading medical device company that delivers proprietary clinical solutions for diagnosing and treating vascular conditions. The company's rapidly growing product line consists of innovative devices across established and emerging areas of coronary and peripheral vascular medicine. For more information, visit [www.vascularsolutions.com](http://www.vascularsolutions.com) and [www.treatveins.com](http://www.treatveins.com).